

Specifications

Visit our website at easystand.com for complete specifications on all EasyStand products and options.

	Evolv medium	Evolv large	Evolv xt	Glider medium	Glider large	StrapStand
Approximate Height Range	4'-5'6" (122-168cm)	5'-6'2" (152-188cm)	6'-6'10" (183-211cm)	4'-5'6" (122-168cm)	5'-6'2" (152-188cm)	5'-6'5" (152-195cm)
Maximum User Weight	200 lbs (90kg)	280 lbs (127kg)	350 lbs (159kg)	200 lbs (90kg)	280 lbs (127kg)	350 lbs (159kg)
Seat Depth Range w/Back (A)	14"-19" (36-48cm)	18"-23" (46-59cm)	20"-25" (50-64cm)	20"-25" (50-64cm)*	24"-29" (61-74cm)*	N/A
Seat to Foot Plate Range (B)	11"-18" (28-46cm)	14"-21" (36-53cm)	17"-23" (43-58cm)	13.5"-16.5" (34-42cm)	15.5"-18.5" (39-47cm)	N/A
Seat Height from Floor	21.5" (55cm)	21.5" (55cm)	23.5" (60cm)	21.5" (55cm)	21.5" (55cm)	N/A
Knee Pad Size	17"Wx9"H (43x23cm)	17"Wx9"H (43x23cm)	17"Wx9"H (43x23cm)	7"Wx7"H (18x18cm)	7"Wx7"H (18x18cm)	17"Wx9"H (43x23cm)
Knee Pad Depth Range	2"-7" (5-18cm)	2"-7" (5-18cm)	2.5"-8" (6-20cm)	N/A	N/A	6"-9.5" (15-24cm)
Weight of Basic Unit	101 lbs (46kg)	104 lbs (47kg)	112 lbs (51kg)	175 lbs (79kg)	178 lbs (81kg)	125 lbs (58kg)
Frame Footprint	26.5"x36" (67x91cm)	26.5"x36" (67x91cm)	28.5"x38" (72x97cm)	26.5"x41" (67x104cm)	26.5"x41" (67x104cm)	28"x34" (71x86cm)